

# DR Shaver Blades and Burs for the CONMED Linvatec Handpiece



## **Our DR SHAVER BLADES AND BURS distinguish themselves with:**

- Very good cutting qualities
- Well seated connections
- Short operative times
- Varying lengths: 85, 130, 155, and 174 mm
- Availability in straight line or 15° angled configurations
- Flexible Blades
- Single use and reusable Blades and Burs

### **DR WAVE CUTTER®**

The DR WAVE CUTTER®, with its double Blade, can effectively and efficiently remove large quantities of tissue.

### **DR PIT BULL®**

The DR PIT BULL®, with three double Blades, is the allrounder for bone and soft tissue. For bursa work and the removal of cruciate ligament remnants. In many instances the Shaver Bur and the Acromionizer need not be used. The attractive solution for clean, rapid, and effective work. Two in one, and you save time and money.

### **DR SHAVER BUR CERATITAN®**

The head and internal channel of the DR BARREL BUR® and the DR ROUND BUR® are clad with Ceratitan (an antibacterial ceramic coating). This coating provides you with superior cut quality as well as very good cutting performance and high contrast for visualization. Under usual circumstances the two aspiration windows prevent “snowstorming”.

**All our DR Shaver Blades and Burs are manufactured at our parent facility in Switzerland – 100% Swiss Made .**

# DR SHAVER BLADES & BURS

Description	Ø (mm)	L (mm)	Ref.No.	Hub
<b>Plastic blue hub</b>	Single use instruments into sterile packaging with 5 pieces inside			
DR WAVE CUTTER	Ø 2.0	85	027.801.230	DROR
DR WAVE CUTTER	Ø 3.0	85	027.301.231	DROR
DR WAVE CUTTER	Ø 3.0	130	027.301.230	DROR
DR WAVE CUTTER	Ø 4.0	130	027.101.230	DROR
DR WAVE CUTTER	Ø 5.0	130	027.201.230	DROR
DR PIT BULL	Ø 3.0	130	027.301.240	DROR
DR PIT BULL	Ø 4.0	130	027.101.240	DROR
DR PIT BULL	Ø 5.0	130	027.201.240	DROR
DR PIT BULL	Ø 7.0	130	027.701.240	DROR
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	027.101.236	DROR
DR BENDING PIT BULL	Ø 3.0	130	027.301.340	DROR
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	027.101.241	DROR
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	027.101.242	DROR
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	027.101.250	DROR
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	027.101.251	DROR
DR BARREL BUR	Ø 4.0	130	137.100.220	DROR
DR BARREL BUR	Ø 5.5	130	137.200.220	DROR
DR ROUND BUR	Ø 4.0	130	127.100.220	DROR
DR ROUND BUR	Ø 5.5	130	127.200.220	DROR
<b>Metal hub silver</b>	Single use instruments into sterile packaging with 5 pieces inside			
DR WAVE CUTTER	Ø 2.0	85	027.801.030	DROR
DR WAVE CUTTER	Ø 3.0	85	027.301.030	DROR
DR WAVE CUTTER	Ø 3.0	130	027.301.031	DROR
DR WAVE CUTTER	Ø 4.0	130	027.101.030	DROR
DR WAVE CUTTER	Ø 5.0	130	027.201.030	DROR
DR PIT BULL	Ø 3.0	85	027.301.047	DROR
DR PIT BULL	Ø 3.0	130	027.301.040	DROR
DR PIT BULL	Ø 4.0	130	027.101.040	DROR
DR PIT BULL	Ø 5.0	130	027.201.040	DROR
DR PIT BULL	Ø 7.0	130	027.701.040	DROR
DR PIT BULL, CONCAVE 15°	Ø 4.0	130	027.101.036	DROR
DR PIT BULL, CONCAVE 15°	Ø 5.0	155	027.201.043	DROR
DR PIT BULL	Ø 4.0	155	027.101.037	DROR
DR PIT BULL	Ø 5.0	155	027.201.037	DROR
DR BENDING PIT BULL	Ø 3.0	130	027.301.140	DROR
DR BENDING TOOL			020.101.103	
DR SYNOVIA PIT BULL	Ø 4.0	130	027.101.041	DROR
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	130	027.101.042	DROR
DR SYNOVIA PIT BULL, CONCAVE 15°	Ø 4.0	174	027.101.050	DROR
DR PIT BULL, CONCAVE 15°	Ø 4.0	174	027.101.051	DROR
DR BARREL BUR	Ø 4.0	130	137.100.020	DROR
DR BARREL BUR	Ø 5.5	130	137.200.020	DROR
DR BARREL BUR	Ø 4.0	174	137.100.021	DROR
DR ROUND BUR	Ø 4.0	130	127.100.020	DROR
DR ROUND BUR	Ø 5.5	130	127.200.020	DROR
DR ROUND BUR	Ø 4.0	174	127.100.021	DROR
DR ROUND BUR	Ø 5.5	174	127.200.021	DROR
<b>Metal hub with gold surface</b>	Reusable instruments into sterile packaging with 1 piece inside			
DR WAVE CUTTER	Ø 4.0	130	027.101.528	DROR
DR WAVE CUTTER	Ø 5.0	130	027.201.528	DROR
DR PIT BULL	Ø 4.0	130	027.101.538	DROR
DR PIT BULL	Ø 5.0	130	027.201.538	DROR
DR PIT BULL	Ø 7.0	130	027.701.538	DROR
DR PIT BULL	Ø 4.0	155	027.101.535	DROR
DR PIT BULL	Ø 5.0	155	027.201.535	DROR
DR SYNOVIA PIT BULL	Ø 4.0	130	027.101.534	DROR
DR BARREL BUR	Ø 4.0	130	137.100.519	DROR
DR BARREL BUR	Ø 5.5	130	137.200.519	DROR
DR ROUND BUR	Ø 4.0	130	127.100.519	DROR
DR ROUND BUR	Ø 5.5	130	127.200.519	DROR

Gold: Can be used 10 times




**DR WAVE CUTTER®**  
(one double blade)




**DR PIT BULL®**  
(three double blades, large window)



**DR SYNOVIA PIT BULL®**  
(three double blades, small window)



**DR BARREL BUR®**



**DR ROUND BUR®**

Our DROR DR SHAVER BLADES & BURS are intended for use with arthroscopy shaver systems and / or handpieces manufactured by CONMED. The firm Conmed Corporation Inc. and DR MEDICAL AG are independent companies.

Version001\_EN